



UNIVERSITY OF CAMBRIDGE INTERNATIONAL EXAMINATIONS  
 General Certificate of Education  
 Advanced Subsidiary Level and Advanced Level

CANDIDATE NAME

CENTRE NUMBER

CANDIDATE NUMBER



**CHEMISTRY** **9701/35**  
 Advanced Practical Skills 1 **October/November 2011**  
**2 hours**

Candidates answer on the Question Paper.  
 Additional Materials: As listed in the Confidential Instructions

**READ THESE INSTRUCTIONS FIRST**

Write your Centre number, candidate number and name on all the work you hand in.  
 Give details of the practical session and laboratory where appropriate, in the boxes provided.  
 Write in dark blue or black pen.  
 You may use a soft pencil for any diagrams, graphs or rough working.  
 Do not use staples, paper clips, highlighters, glue or correction fluid.  
 DO **NOT** WRITE IN ANY BARCODES.

Answer **all** questions.  
 You may lose marks if you do not show your working or if you do not use appropriate units.  
 Use of a Data Booklet is unnecessary.

Qualitative Analysis Notes are printed on pages 11 and 12.


At the end of the examination, fasten all your work securely together.  
 The number of marks is given in brackets [ ] at the end of each question or part question.

<b>Session</b>
<b>Laboratory</b>

For Examiner's Use	
1	
2	
<b>Total</b>	

This document consists of **10** printed pages and **2** blank pages.

- 1 You are required to investigate the effect of temperature on the rate of reaction of peroxydisulfate ions with iodide ions. Iodide ions are oxidised to iodine by peroxydisulfate ions.



**FA 1** is aqueous potassium peroxydisulfate,  $\text{K}_2\text{S}_2\text{O}_8$ .

**FA 2** is an aqueous solution containing a mixture of potassium iodide, KI, sodium thiosulfate,  $\text{Na}_2\text{S}_2\text{O}_3$ , and starch.

When **FA 1** and **FA 2** are mixed together the potassium peroxydisulfate reacts with the potassium iodide to make iodine. As soon as this iodine is formed, it reacts with the sodium thiosulfate and is turned back into iodide ions. Only when all the sodium thiosulfate has reacted does iodine remain in the solution. The solution then turns blue-black because of the presence of the starch indicator. The rate of reaction can be determined by the time it takes for a blue-black colour to first appear in the colourless mixture.

### (a) Method


**Read through the method and prepare a table on page 3 to record the initial and final temperatures and the reaction time for each experiment, before starting any practical work.**

- Half-fill a  $250\text{cm}^3$  beaker with water to act as a water bath.
- Place it on a tripod and gauze and heat it with a Bunsen burner to about  $65^\circ\text{C}$  then remove the Bunsen. While your water is being heated continue with the following steps of the method.
- Fill the burette, labelled **FA 1**, with the aqueous potassium peroxydisulfate, **FA 1**.
- Fill the burette, labelled **FA 2**, with the mixture of solutions, **FA 2**.
- Measure  $10.0\text{cm}^3$  of **FA 1** into a boiling tube.
- Measure  $10.0\text{cm}^3$  of **FA 2** into a second boiling tube.
- Place both boiling tubes in the water bath.
- Clamp one of the tubes and place a thermometer in this tube.
- When the temperature of this solution has reached about  $60^\circ\text{C}$ , pour the contents of the second tube into the clamped tube. Start timing immediately, note the temperature and stir the mixture.
- Record this initial temperature.
- Stop timing as soon as the blue-black colour appears. Record this reaction time to the **nearest second** and record the final temperature.
- Repeat the experiment at decreasing temperatures as many times as necessary to generate data for plotting a graph. The experiment should not be performed at a temperature below about  $30^\circ\text{C}$ . The temperature of the water bath may be adjusted by adding cold water or by reheating. (Boiling tubes may be rinsed and reused.)



- (c) (i) Using your values in (b), plot a graph of 'rate' ( $y$ -axis) against mean temperature ( $x$ -axis). Choose suitable scales to allow you to extrapolate the graph to include the 'rate' at  $20^{\circ}\text{C}$ .

For  
Examiner's  
Use



I	
II	
III	
IV	

- (ii) Use your graph to obtain a value for the 'rate' at  $20^{\circ}\text{C}$ .

'rate' at  $20^{\circ}\text{C}$  is ..... [4]

(d) It has been suggested that an increase in temperature of 10°C will double the rate of reaction. Use two pairs of temperatures from your graph to confirm or deny this statement.

.....  
.....  
.....  
..... [2]

(e) (i) Which of your experiments has the greatest percentage error in timing?

.....

(ii) Calculate the percentage error in (i).  
You may assume that the error in measuring the time for a reaction is  $\pm 0.5$  seconds.

error in (i) = ..... % [2]

(f) A student had difficulty in drawing a line of best fit. Identify a source of error in the experimental **procedure**. Do not include any errors involving the precision of apparatus.

.....  
..... [1]

(g) Suggest a modification that could be used to reduce this error.

.....  
..... [1]

(h) Using **FA 1**, **FA 2** and distilled water, describe how you could investigate the effect of **concentration** of potassium peroxydisulfate on the rate of reaction.

.....  
.....  
.....  
.....  
..... [3]

[Total: 24]

## 2 Qualitative analysis

For  
Examiner's  
Use

At each stage of any test you are to record details of the following.

- colour changes seen
- the formation of any precipitate
- the solubility of such precipitates in an excess of the reagent added

Where gases are released they should be identified by a test, **described in the appropriate place in your observations.**

You should indicate clearly at what stage in a test a change occurs.  
Marks are **not** given for chemical equations.

**No additional tests for ions present should be attempted.**

**If any solution is warmed, a boiling tube MUST be used.**

**Rinse and reuse test-tubes and boiling tubes where possible.**

**Where reagents are selected for use in a test the full name or correct formula of the reagents must be given.**

**FA 3, FA 4 and FA 5** are aqueous solutions containing one cation and one anion. One of these solutions is a dilute acid and this is the only acid present.

- (a) (i) Select a single chemical reagent from those supplied which would allow you to identify the dilute acid. You may not use indicator paper.

reagent .....

- (ii) Use this reagent to test all three solutions and record your results in an appropriate form in the space below.

I	
II	
III	
IV	
V	

- (iii) From your observations in (ii), identify which solution is the dilute acid.

**FA** ..... is the dilute acid.

[5]

- (b) The acid you have identified in (a)(iii) is dilute sulfuric acid. Complete the following table.

For  
Examiner's  
Use

<i>test</i>	<i>observations</i>
To 1 cm depth of <b>FA 3</b> in a test-tube, add 1 cm depth of <b>FA 4</b> , then	
add excess hydrochloric acid.	
To 1 cm depth of <b>FA 4</b> in a test-tube, add 1 cm depth of <b>FA 5</b> , then	
add excess hydrochloric acid.	
To 1 cm depth of <b>FA 5</b> in a test-tube, add 1 cm depth of <b>FA 3</b> , then	
add excess hydrochloric acid.	

I	
II	
III	
IV	

[4]

- (c) For the two unidentified solutions, complete the following table.

<i>test</i>	<i>observations</i>	
	<b>FA .....</b>	<b>FA .....</b>
To 1 cm depth of unknown in a boiling tube, add NaOH(aq)		
warm the tube carefully		

[2]

(d) From your observations in (a), (b) and (c), identify the ions present in the two solutions tested in (c), giving the relevant evidence for each. If you have not been able to identify one or more of the ions, explain why the evidence obtained was insufficient.

For  
Examiner's  
Use

FA ..... cation ..... evidence .....

.....

.....

anion ..... evidence .....

.....

.....

FA ..... cation ..... evidence .....

.....

.....

anion ..... evidence .....

.....

.....

[4]

(e) If one of the aqueous anions was a bromide, what would be the minimum evidence needed for its identification?

.....

.....

..... [1]

[Total: 16]

I	
II	
III	
IV	



**BLANK PAGE**

**BLANK PAGE**

---

Permission to reproduce items where third-party owned material protected by copyright is included has been sought and cleared where possible. Every reasonable effort has been made by the publisher (UCLES) to trace copyright holders, but if any items requiring clearance have unwittingly been included, the publisher will be pleased to make amends at the earliest possible opportunity.

University of Cambridge International Examinations is part of the Cambridge Assessment Group. Cambridge Assessment is the brand name of University of Cambridge Local Examinations Syndicate (UCLES), which is itself a department of the University of Cambridge.

## Qualitative Analysis Notes

Key: [ppt. = precipitate]

## 1 Reactions of aqueous cations

ion	reaction with	
	NaOH(aq)	NH <sub>3</sub> (aq)
aluminium, Al <sup>3+</sup> (aq)	white ppt. soluble in excess	white ppt. insoluble in excess
ammonium, NH <sub>4</sub> <sup>+</sup> (aq)	no ppt. ammonia produced on heating	–
barium, Ba <sup>2+</sup> (aq)	no ppt. (if reagents are pure)	no ppt.
calcium, Ca <sup>2+</sup> (aq)	white ppt. with high [Ca <sup>2+</sup> (aq)]	no ppt.
chromium(III), Cr <sup>3+</sup> (aq)	grey-green ppt. soluble in excess giving dark green solution	grey-green ppt. insoluble in excess
copper(II), Cu <sup>2+</sup> (aq)	pale blue ppt. insoluble in excess	blue ppt. soluble in excess giving dark blue solution
iron(II), Fe <sup>2+</sup> (aq)	green ppt. turning brown on contact with air insoluble in excess	green ppt. turning brown on contact with air insoluble in excess
iron(III), Fe <sup>3+</sup> (aq)	red-brown ppt. insoluble in excess	red-brown ppt. insoluble in excess
lead(II), Pb <sup>2+</sup> (aq)	white ppt. soluble in excess	white ppt. insoluble in excess
magnesium, Mg <sup>2+</sup> (aq)	white ppt. insoluble in excess	white ppt. insoluble in excess
manganese(II), Mn <sup>2+</sup> (aq)	off-white ppt. rapidly turning brown on contact with air insoluble in excess	off-white ppt. rapidly turning brown on contact with air insoluble in excess
zinc, Zn <sup>2+</sup> (aq)	white ppt. soluble in excess	white ppt. soluble in excess

[Lead(II) ions can be distinguished from aluminium ions by the insolubility of lead(II) chloride.]

## 2 Reactions of anions

<i>ion</i>	<i>reaction</i>
carbonate, $\text{CO}_3^{2-}$	$\text{CO}_2$ liberated by dilute acids
chromate(VI), $\text{CrO}_4^{2-}(\text{aq})$	yellow solution turns orange with $\text{H}^+(\text{aq})$ ; gives yellow ppt. with $\text{Ba}^{2+}(\text{aq})$ ; gives bright yellow ppt. with $\text{Pb}^{2+}(\text{aq})$
chloride, $\text{Cl}^-(\text{aq})$	gives white ppt. with $\text{Ag}^+(\text{aq})$ (soluble in $\text{NH}_3(\text{aq})$ ); gives white ppt. with $\text{Pb}^{2+}(\text{aq})$
bromide, $\text{Br}^-(\text{aq})$	gives cream ppt. with $\text{Ag}^+(\text{aq})$ (partially soluble in $\text{NH}_3(\text{aq})$ ); gives white ppt. with $\text{Pb}^{2+}(\text{aq})$
iodide, $\text{I}^-(\text{aq})$	gives yellow ppt. with $\text{Ag}^+(\text{aq})$ (insoluble in $\text{NH}_3(\text{aq})$ ); gives yellow ppt. with $\text{Pb}^{2+}(\text{aq})$
nitrate, $\text{NO}_3^-(\text{aq})$	$\text{NH}_3$ liberated on heating with $\text{OH}^-(\text{aq})$ and <i>Al</i> foil
nitrite, $\text{NO}_2^-(\text{aq})$	$\text{NH}_3$ liberated on heating with $\text{OH}^-(\text{aq})$ and <i>Al</i> foil, $\text{NO}$ liberated by dilute acids (colourless $\text{NO} \rightarrow$ (pale) brown $\text{NO}_2$ in air)
sulfate, $\text{SO}_4^{2-}(\text{aq})$	gives white ppt. with $\text{Ba}^{2+}(\text{aq})$ or with $\text{Pb}^{2+}(\text{aq})$ (insoluble in excess dilute strong acids)
sulfite, $\text{SO}_3^{2-}(\text{aq})$	$\text{SO}_2$ liberated with dilute acids; gives white ppt. with $\text{Ba}^{2+}(\text{aq})$ (soluble in excess dilute strong acids)

## 3 Tests for gases

<i>gas</i>	<i>test and test result</i>
ammonia, $\text{NH}_3$	turns damp red litmus paper blue
carbon dioxide, $\text{CO}_2$	gives a white ppt. with limewater (ppt. dissolves with excess $\text{CO}_2$ )
chlorine, $\text{Cl}_2$	bleaches damp litmus paper
hydrogen, $\text{H}_2$	“pops” with a lighted splint
oxygen, $\text{O}_2$	relights a glowing splint
sulfur dioxide, $\text{SO}_2$	turns potassium dichromate(VI) (aq) from orange to green